

TabTime Vibrating Watch

Note there are **two buttons** to adjust the settings on this watch. As you place the watch on your wrist, locate these two button by feeling the side of the watch face. The buttons are tiny and situated about a centimetre apart on the same side (thin edge of the watch face).

Setting the Time

Tap the **right button** to display the current time.

Now hold the **right button** down for 3 seconds.

The **left** button **increases** the hours and the **right** button **decreases** the hours.

Hold the right button for 3 seconds **to confirm** the hour.

Now you can set the minutes.

Hold the right button for 3 seconds **to confirm** the minutes.

Activating the Alarm

Use the **left button** to cycle through the 10 alarms.

Press the **right button** to turn on/off the alarm.

Changing the Alarm Time

Tap the **left button** to display the alarm time, then hold the left button, displayed time should flash.

Now set the hour for your desired alarm time.

Hold the **left button** for 3 seconds to confirm the hour.

Now you can set the minutes.

Hold the **left button** for 3 seconds to confirm the minutes.

Discrete Vibrating Alarms

Alarm “A0” vibrates for 40 seconds.

Alarm “A1” to “A9” vibrate for 20 seconds.

Press the right button to stop the Alarm.

Charging the Watch

The strap is easy to disconnect and reconnect.

Revealing the charging port, small USB charging cable included.

Only takes two hours to charge fully.